THE STEP TECHNOLOGY (Space & Time Extinction Profile)

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The measurement cycle to extrapolate the shelf life for the product at room temperature during three years, is shorter and take only 15 hours.

The principle (Fig.1) is based on the transmission of a light signal close to infra red, coupled with a centrifugation system with a centrifugal force 1000g at a rotation speed of 3000rpm at a temperature of 35°C (1).

Software (SEPView) records the transmission intensity of a sample in terms of time and position in order to draw a "transmission profile" (Fig.2).

The additional advantages from this instrumental method, apart from being time-saving, are: the small quantity needed to analyze: only 0,5ml, and the possibility to assess eight samples during one cycle.

The profile interpretation is immediately obtained and the creation of a references data base is not needed. The results are immediate, accurate and could be correlated by visual observation of measurement cells.

Neither the "benchmark" nor the "prototype" showed changes in appearance, in color, in odor or in texture.

Thus, the products seem stable.

Texture assessment:

We have also two kinds of test to characterize formula textures:

a subjective test based on panel answers (2), and an objective test based on measurement by a Texture Analyzer TA-XT Plus 5kg (Stable MicroSystem) (3).

Sensory assessment is ore and more use by cosmetic companies to characterize emulsions and it would seem to be a relevant tool to complete the development procedure for new skin care products.

Each formula is proposed to a panel of ten women between 25-35 years old to assess their sensory properties, I.e. their respective texture and skin feel (Coty procedure).

Esquema de funcionamiento



Fig. 2 The classic transmission profile

Reference sample	Sedimentation Calculation		Extrapolation 3 Years	Observation cell
	SLOPE (µm/year)	SLOPE (mm/year)		
Prototype	0,0044	0,1368	medium sedimentation, not visible	ok
Benchmark	0,0040	0,1244	medium sedimentation, not visible	ok

Table 2: The LumFuge data

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Stability Analyser
LUMiFuge®

- ✓Very slow separation & demixing pr. (month-years)
- ✓Very stable, very high viscous, very high concentrations
- ✓High throughput with 8 samples measured at once
- ✓SHELF LIFE, STABILITY

Dispersion Analyser
LUMiSizer®

- ✓Very slow separation & demixing pr. (month-years)
- ✓Very stable, very high viscous, very high concentrations
- ✓High throughput with 12 samples measured at once
- ✓PSD, Very small particles or droplets (0.05 µm – 1000µm)
- ✓SHELF LIFE, STABILITY, PSD

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

Separation Analyser
LUMiReader®

- ✓Fast separation / demixing processes (mins – hours)
- ✓Fragile / shear sensitive samples
- ✓Analysis under original conditions
- ✓Separation enhancement by inclination up to 10g
- ✓COALESCENCE, FLOCCULATION

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Eye-Nose-Cream development	Sedimentation velocity at 2300 xg [µm/sec]	Sedimentation velocity at gravity [µm/d]
A (blue line)	0.109	4.1
B (light green line)	0.111	4.2
Ref. market (red line)	0.047	1.8
C (dark green line)	0.033	1.3

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